

- 8 Read the passage and answer the questions that follow.

X-Rays

A modern X-ray tube is illustrated in Fig. 8.1. A heated cathode produces electrons, known as cathode ray, and these are accelerated by the potential difference between the cathode and anode. X-ray tubes used by dentists typically use potential differences of 60 kV. As the electrons hit the anode, X-rays are produced. Tungsten metal is an ideal material for the target, or anode, of an X-ray machine. In practice, only about 1 % of the energy of the electron beam produces X-rays. The rest is wasted as heat.




Fig. 8.1

The X-ray spectrum is both a continuous and a superimposed line spectrum, as shown in Fig. 8.2. Fig. 8.2a shows the X-ray spectrum and Fig. 8.2b shows how X-ray intensity varies with the wavelength. The wavelengths of the X-ray spectra are determined by the element from which the target is made.




Fig. 8.2

Electrons in the tungsten atom have energy levels in a similar way to the energy levels for a hydrogen atom. Fig. 8.3 shows the logarithm of the energy E of some of the energy levels of electrons in the tungsten atom. E is measured in electron-volts (eV).




Fig. 8.3

Attenuation of X-ray

When X-ray passes through a medium, it is absorbed by the medium.

The intensity I of a beam of X-ray is related to the distance x it travels through a medium by the relationship:

$$I = I_0 e^{-\mu x}$$

where I_0 is the incident intensity and μ is the **total linear attenuation coefficient**, a constant which depends on the medium and the photon energy of the X-rays.

The decrease variation with distance of the intensity is shown in Fig. 8.4.




Fig. 8.4

The penetrating power, or quality, of a radiation can conveniently be described in terms of the thickness of material needed to reduce the intensity to half the original value. This is called the **half-value thickness (HVT)**, $x_{\frac{1}{2}}$.

- (a) In the production of X-rays, electrons are emitted from the cathode and accelerated in an electric field. The electron energy is often described in terms of electron-volts.

Explain what is meant by the term *electron-volt*.

.....

 [2]

- (b) (i) When an electron strikes the target anode, it can emit energy in creating a single photon.

Determine the minimum wavelength for the X-ray spectrum used by dentist.

minimum wavelength = m [3]

- (ii) Suggest how the continuous spectrum of X-rays in Fig. 8.2 are produced.

.....

 [1]

- (c) An electron falling from the L to the K level gives rise to the K_α line. The photon energy of electrons falling from level L to level K depends on the element used for the target. The table in Fig. 8.5 shows the photon energy and the proton number Z for three elements.

element	proton number Z	K_α photon energy / keV
chromium	24	5.40
copper	29	8.03
silver	47	22.10
tungsten	74	

Fig. 8.5

In 1914 Henry G.J. Moseley discovered empirically that the wavelength λ of the K_α line is related to the proton number by the relationship:

$$\sqrt{\frac{1}{\lambda}} \propto Z$$

which is known as the Moseley's law.

- (i) Determine the K_α photon energy for the tungsten atom.

energy = keV [2]

- (ii) Determine quantitatively, without drawing a graph, whether Moseley's law is valid.

[3]

- (iii) Explain why the K_β line has a shorter wavelength than the K_α line.

.....

 [1]

- (d) Suggest and explain one feature of tungsten that makes it an ideal metal for the target of an X-ray machine.

.....

 [2]

- (e) A beam of X-rays is incident on a hand and the X-ray photograph is taken, as shown in Fig. 8.6.



Fig. 8.6

- (i) The beam is attenuated as it passes through the bone. The total linear attenuation coefficient of the bone is 0.528 cm^{-1} .

1. Show that the penetrating power of the X-ray in human bone is 1.31 cm

[1]

2. Calculate the depth of bone required to reduce the intensity of the beam to 40.0% of its initial value.

depth of bone = cm [2]

- (ii) An identical beam of X-rays is also attenuated by 40% as it passes through 3.87 cm of human muscle.
 Calculate the total linear attenuation coefficient of human muscle.

total linear attenuation coefficient = cm^{-1} [1]

[Turn over

- (iii) Briefly explain how these different values of attenuation coefficient help in the interpretation of X-ray pictures taken of human bodies.

.....

.....

..... [2]